
Product Data Sheet

Product Name: H-Lys-Lys-Lys-OH

Cat. No.: GA22698

Chemical Properties

Cas. No. 13184-14-0

Formula $C_{18}H_{38}N_6O_4$ M.Wt 402.54

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Trilysine has found use in medicine as component of polymeric hydrogel adhesives for wound closure and other surgical applications. Such hydrogels form in situ and can also serve as depot for the slow release of drugs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA