
Product Data Sheet

Product Name: H-Met-Gly-Met-Met-OH

Cat. No.: GA22731

Chemical Properties

Cas. No. 14517-45-4

Formula $C_{17}H_{32}N_4O_5S_3$ M.Wt 468.66

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for human methionine aminopeptidase 2, Km 0.34 ± 0.04 mM and kcat 170 ± 6 min⁻¹ at pH 7.5 and 30°C.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
