
Product Data Sheet

Product Name: HOBT-6-sulfonamidomethyl resin . HCl (200-400 mesh, 0.6-0.9 mmol/g)

Cat. No.: GA22760

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Highly efficient polymeric reagent for the synthesis of amides. Simple filtration allows the separation of the product from the polymer. Can be used in high throughput screening assays. Devocelle et al. employed it in the parallel synthesis of hydroxam

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA