
Product Data Sheet

Product Name: H-Phe-Arg-OH

Cat. No.: GA22788

Chemical Properties

Cas. No. 1238-09-1

Formula $C_{15}H_{23}N_5O_3$ M.Wt 321.38

Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$, protect from light

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Phe-Arg is cleaved from bradykinin by angiotensin-converting enzyme (ACE). See also H-Ser-Pro-OH and Bradykinin (1-5) for further cleavage products.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA