

---

**Product Data Sheet**

---

Product Name: H-Phe-Gly-Gly-OH

Cat. No.: GA22793

**Chemical Properties**

Cas. No. 23576-42-3

Formula  $C_{13}H_{17}N_3O_4$  M.Wt 279.3Solubility  $H_2O$  : 25 mg/mL (89.51 mM; Need ultrasonic); DMSO : 2 mg/mL (7.16 mM; ultrasonic and warming and heat to 60°C) Store Storage at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Substrate for the characterization of tripeptide aminopeptidase from monkey brain.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA