
Product Data Sheet

Product Name: H-Phe-Pro-Arg-OH

Cat. No.: GA22811

Chemical Properties

Cas. No. 37553-80-3

Formula $C_{20}H_{30}N_6O_4$ M.Wt 418.5

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

FPR is a tripeptide analog of the thrombin-receptor-derived peptide SFLLR. In the chick chorioallantoic membrane system it inhibited angiogenesis and caused a complete reversal of the angiogenesis-promoting effect of the thrombin-receptor-activating te

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA