
Product Data Sheet

Product Name: H-Pro-Thr-His-Ile-Lys-Trp-Gly-Asp-OH

Cat. No.: GA22870

Chemical Properties

Cas. No. 117620-76-5

Formula C₄₄H₆₄N₁₂O₁₂ M.Wt 953.07

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Tuna Muscle Peptide, PTHIKWGD, an angiotensin I-converting enzyme (ACE) inhibitor originally isolated from tuna muscle.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA