
Product Data Sheet

Product Name: H-Ser-Gln-Asn-Tyr-Pro-Ile-Val-OH

Cat. No.: GA22886

Chemical Properties

Cas. No. 118071-31-1

Formula C₃₇H₅₇N₉O₁₂ M.Wt 819.91

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The heptapeptide SQNYPIV has a sequence similar to the known p17-p24 junction in the gag polyprotein and is efficiently used as a substrate by the HIV-1 protease. SQNYPIV binds to the tetanus C fragment (TetC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA