
Product Data Sheet

Product Name: H-Tyr(SO₃H)-OH

Cat. No.: GA22961

Chemical Properties

Cas. No. 956-46-7

Formula C₉H₁₁NO₆S M.Wt 261.26

Solubility Soluble in DMSO Storage Store at -15°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sulfotyrosine has been used by Liu et al. for conveniently obtaining tyrosine-sulfated proteins by expression in E. coli with an expanded genetic code.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
