
Product Data Sheet

Product Name: H-Tyr-Cys-Phe-Ala-Trp-Lys-Thr-Phe-Cys-OH

Cat. No.: GA22966

Chemical Properties

Cas. No. 647016-22-6

Formula C₅₇H₇₁N₁₁O₁₂S₂

M.Wt

1166.39

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Selective somatostatin receptor subtype sst4 agonist. In competitive radioligand displacement experiments using the non-selective radioligand ([¹²⁵I]-(Leu,D-Trp²²,Tyr²)-Somatostatin-28, YCFAWKTFC exhibited a high affinity for ss

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA