
Product Data Sheet

Product Name: H-Val-His-Leu-Thr-Pro-OH

Cat. No.: GA23010

Chemical Properties

Cas. No. 93913-38-3

Formula C₂₆H₄₃N₇O₇ M.Wt 565.67

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

N-terminal pentapeptide of human hemoglobin beta (Hb A and and the pathological HbS, which differ in position 6). Cellulose membrane-bound VHLTP was used as ligand for the anti-HbA1c antibody.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA