

---

**Product Data Sheet**

---

Product Name: Hymenistatin I

Cat. No.: GA23037

**Chemical Properties**

Cas. No. 129536-23-8

Formula C<sub>47</sub>H<sub>72</sub>N<sub>8</sub>O<sub>9</sub> M.Wt 893.14

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Hymenistatin I (HS-I) exerts an immunosuppressive effect in the humoral and cellular immune responses comparable with that of cyclosporin A (CsA). Comparison of the influence of HS-I and CsA on cytokine production suggests that the mechanisms of the in

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA