
Product Data Sheet

Product Name: IFN- γ Antagonist

Cat. No.: GA23098

Chemical Properties

Cas. No. 158040-83-6

Formula $C_{115}H_{194}N_{34}O_{34}S$

M.Wt 2629.08

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This peptide derived from the ligand binding site of the human γ -interferon (IFN- γ) receptor is an antagonist of human IFN- γ . It inhibited human IFN- γ -induced expression of HLR/DR antigen on Colo 205 cells with an approximate IC50 o

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA