

---

**Product Data Sheet**

---

Product Name: Iodopolystyrene resin (200-400 mesh, &gt; 1 mmol/g)

Cat. No.: GA23109

**Chemical Properties**

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

This "polymeric iodide" has been oxidized to yield iodosyl compounds such as the iodosyl acetate (M.L.Hallensleben). Being an excellent leaving group, iodine may also be replaced with a metal (e.g. lithium). A broad range of resin derivatives

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA