
Product Data Sheet

Product Name: Isovaleryl-Phe-Nle-Sta-Ala-Sta-OH

Cat. No.: GA23110

Chemical Properties

Cas. No. 98092-14-9

Formula C₃₉H₆₅N₅O₉ M.Wt 747.97

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The pepstatin analog SR 42128 Inhibits renin activity.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA