
Product Data Sheet

Product Name: JKC-301
Cat. No.: GA23112

Chemical Properties

Cas. No. 136553-96-3

Formula $C_{32}H_{44}N_6O_7$ M.Wt 624.74

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

JKC-301, a highly potent ET-A receptor antagonist.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA