
Product Data Sheet

Product Name: Lepirudin

Cat. No.: GA23131

Chemical Properties

Cas. No. 138068-37-8

Formula $C_{287}H_{440}N_{80}O_{111}S_6$

M.Wt 6979.52

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Synthetic peptide. Lepirudin is an anticoagulant which functions as a direct thrombin inhibitor.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA