
Product Data Sheet

Product Name: Leptin (116-130) (human)

Cat. No.: GA23132

Chemical Properties

Cas. No. 2243207-12-5

Formula C₇₀H₁₀₆N₁₈O₂₄S

M.Wt 1615.76

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Fragment of human leptin, an adipose tissue-derived protein encoded by the obese gene. Leptin has been implicated in the control of food intake, energy expenditure and reproduction. Additional roles in the regulation of immunity, inflammation, and hema

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA