
Product Data Sheet

Product Name: LL-37 FKR

Cat. No.: GA23145

Chemical Properties

Cas. No. 913736-94-4

Formula $C_{118}H_{199}N_{37}O_{30}$

M.Wt 2616.11

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

LL-37 FKR (FK-21), C-terminal fragment of LL-37, showed antimicrobial activity and similar concentration-dependent chemotactic activities on granulocytes as the full-length peptide.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA