
Product Data Sheet

Product Name: LL-37 KRI

Cat. No.: GA23147

Chemical Properties

Cas. No. 672333-78-7

Formula C₁₀₉H₁₉₀N₃₆O₂₉

M.Wt

2468.93

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

KR-20 corresponds to LL-37 (18-37). The fragment, which has been detected in human sweat, shows increased antimicrobial activity.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA