

---

**Product Data Sheet**

---

Product Name: LL-37 LLG

Cat. No.: GA23148

**Chemical Properties**

Cas. No. 740800-39-9

Formula  $C_{119}H_{196}N_{32}O_{29}$ 

M.Wt 2539.06

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA