

---

**Product Data Sheet**

---

Product Name: LL-37 RKS

Cat. No.: GA23149

**Chemical Properties**

Cas. No. 672345-61-8

Formula  $C_{169}H_{292}N_{54}O_{45}$ 

M.Wt

3800.51

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The antimicrobially active LL-37 fragments RK-31 and KR-20 have been identified as physiological components of human sweat.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA