
Product Data Sheet

Product Name: LQEQ-19 (mouse, rat)

Cat. No.: GA23152

Chemical Properties

Cas. No. 322644-72-4

Formula $C_{106}H_{165}N_{29}O_{34}$

M.Wt

2389.65

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

LQEQ-19 (VGF599-617) is a Pro-VGF-derived peptide that induces penile erection when injected into the paraventricular nucleus of male rats. The effect is mediated by increased nitric oxide production in the hypothalamic nucleus causing activation of pa

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA