

---

**Product Data Sheet**

---

Product Name: Melanocyte Protein PMEL 17 (44-59) (human, bovine, mouse)

Cat. No.: GA23187

**Chemical Properties**

Cas. No. 195523-86-5

Formula  $C_{95}H_{137}N_{27}O_{28}$  M.Wt 2105.3

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA