
Product Data Sheet

Product Name: Mitogenic Pentapeptide

Cat. No.: GA23202

Chemical Properties

Cas. No. 87173-03-3

Formula $C_{67}H_{124}N_6O_{14}S$ M.Wt 1269.82

Solubility Soluble in DMSO Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

The palmitoylated pentapeptide Pam3Cys-Ser-Ser-Asn-Ala represents a potent activator for monocytes/macrophages and B lymphocytes.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA