
Product Data Sheet

Product Name: Myelin Oligodendrocyte Glycoprotein (35-55) (human)

Cat. No.: GA23207

Chemical Properties

Cas. No. 163158-19-8

Formula $C_{120}H_{179}N_{35}O_{28}S$

M.Wt 2592.03

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

In contrast to the encephalitogenic rat homolog, myelin oligodendrocyte glycoprotein (35-55) (human) does not appear to be a major T-cell epitope in multiple sclerosis probably due to its lower amount of helical structure formed when dissolved in micel

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA