
Product Data Sheet

Product Name: Myelin Oligodendrocyte Glycoprotein (35-55) amide (rat, mouse)
Cat. No.: GA23208

Chemical Properties

Cas. No. 2022956-48-3

Formula $C_{118}H_{178}N_{36}O_{28}S$ M.Wt 2581

Solubility H_2O : ≥ 100 mg/mL (38.75 mM) Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA