
Product Data Sheet

Product Name: Neuropeptide Y (3-36) (porcine)

Cat. No.: GA23258

Chemical Properties

Cas. No. 143863-88-1

Formula $C_{176}H_{271}N_{53}O_{54}$

M.Wt

3993.41

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This molecular form of NPY was found in porcine brain, which in contrast to NPY binds selectively to Y receptors.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA