
Product Data Sheet

Product Name: Octaarginine

Cat. No.: GA23316

Chemical Properties

Cas. No. 148796-86-5

Formula C₄₈H₉₈N₃₂O₉

M.Wt

1267.52

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The cell-penetrating peptide octa-arginine is a potent inhibitor of proteasome activities.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA