
Product Data Sheet

Product Name: Osteoblast-Adhesive Peptide

Cat. No.: GA23325

Chemical Properties

Cas. No. 193613-75-1

Formula $C_{21}H_{43}N_{11}O_6$ M.Wt 545.64

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

KRSR selectively enhanced heparan sulfate-mediated osteoblast adhesion mechanisms. Elucidating such peptide sequence is of importance to the field of bone-cell/tissue engineering as well as to the design of proactive dental/orthopedic biomaterials. Note

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA