
Product Data Sheet

Product Name: Ovalbumin (154-159)

Cat. No.: GA23333

Chemical Properties

Cas. No. 1370698-94-4

Formula $C_{28}H_{52}N_{10}O_9$ M.Wt 672.78

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The hexapeptide TNGIIR exhibited higher inhibition of ACE than several other peptides tested in comparison, as found by Yu et al. The concentration of the peptide, necessary to inhibit 50% of the activity of ACE, was 70 μ M.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA