
Product Data Sheet

Product Name: Ovotransferrin (328-332)

Cat. No.: GA23335

Chemical Properties

Cas. No. 1226776-54-0

Formula $C_{25}H_{46}N_8O_7$ M.Wt 570.69

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

ACE-inhibitory peptide associated with antioxidant and anticoagulation properties, IC50 20 μM.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA