
Product Data Sheet

Product Name: Oxyntomodulin (30-37) (bovine, dog, porcine)

Cat. No.: GA23336

Chemical Properties

Cas. No. 81117-26-2

Formula $C_{39}H_{72}N_{16}O_{12}$ M.Wt 957.1

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

The octapeptide KA8 corresponds to the C-terminal fragment of bovine glucagon-37 (oxyntomodulin). The basic octapeptide was shown to be biologically active.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA