
Product Data Sheet

Product Name: PACAP Related Peptide (1-29) (rat)

Cat. No.: GA23339

Chemical Properties

Cas. No. 132769-35-8

Formula $C_{148}H_{242}N_{42}O_{45}S$

M.Wt

3361.87

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The PACAP precursor contains another peptide called PACAP-related peptide, which is expressed in the parvocellular part of the paraventricular nucleus (PVN), implying its role as hypothalamic regulatory peptide. So far no information on this peptide in

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA