
Product Data Sheet

Product Name: PACAP-38 (28-38) (human, chicken, mouse, ovine, porcine, rat)

Cat. No.: GA23340

Chemical Properties

Cas. No. 160489-86-1

Formula $C_{61}H_{110}N_{24}O_{14}$ M.Wt 1403.7

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA