
Product Data Sheet

Product Name: PAR-2 (1-6) (human)

Cat. No.: GA23348

Chemical Properties

Cas. No. 202933-49-1

Formula $C_{28}H_{53}N_7O_8$ M.Wt 615.77

Solubility DMSO : 25 mg/mL (40.60 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The PAR-2 agonist SLIGKV can be used for the investigation of receptor functions.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA