
Product Data Sheet

Product Name: PAR-3 (1-6) amide (human)

Cat. No.: GA23350

Chemical Properties

Cas. No. 1872435-09-0

Formula C₂₉H₄₆N₁₀O₇ M.Wt 646.75

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

TFRGAP-amide, human PAR-3-derived tethered ligand sequence which does not activate PAR-3 but rather activates PAR-1 and PAR-2, either in Jurkat or in other PAR-expressing cells.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA