
Product Data Sheet

Product Name: Peptide WE-14

Cat. No.: GA23356

Chemical Properties

Cas. No. 115136-18-0

Formula $C_{72}H_{116}N_{18}O_{24}S$

M.Wt

1649.89

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

WE-14, a chromogranin A-derived peptide, was isolated from a human ileal carcinoid tumor. Its sequence WSKMDQLAKELTAE is highly conserved in many mammalian species and flanked by typical processing sites. Such factors would indicate a potential physiolo

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA