
Product Data Sheet

Product Name: PHI-27 (porcine)

Cat. No.: GA23367

Chemical Properties

Cas. No. 80458-29-3

Formula $C_{136}H_{216}N_{36}O_{40}$

M.Wt 2995.43

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Porcine intestinal peptide which belongs to the glucagon-secretin family. It exhibits biological activities similar to vasoactive intestinal peptide (VIP) and secretin.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA