

---

**Product Data Sheet**

---

Product Name: Preangiotensinogen (11-14) (human)

Cat. No.: GA23375

**Chemical Properties**

Cas. No. 102813-98-9

Formula  $C_{21}H_{35}N_7O_6$  M.Wt 481.55

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Cleavage product of human preangiotensinogen (1-14) (M-1120) by the action of renin.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA