
Product Data Sheet

Product Name: Prepro-Atrial Natriuretic Factor (26-55) (human)

Cat. No.: GA23378

Chemical Properties

Cas. No. 112160-82-4

Formula $C_{152}H_{236}N_{38}O_{51}S_3$

M.Wt 3507.97

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Prepro-atrial natriuretic factor (26-55) (human) (NT-pro-ANP (1-30)) is one of four peptide hormones derived from the atrial natriuretic prohormone. It is also known as long acting natriuretic peptide. The peptide has potent vasodilatory properties equ

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA