
Product Data Sheet

Product Name: Prepro-Atrial Natriuretic Factor (56-92) (human)

Cat. No.: GA23379

Chemical Properties

Cas. No. 112199-06-1

Formula $C_{173}H_{270}N_{44}O_{57}$

M.Wt

3878.31

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

NT-pro-ANP (31-67) is used for raising antibodies.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA