

---

## Product Data Sheet

---

Product Name: Prepro-Killer Cell Immunoglobulin-Like Receptor 3DL3 (23-31)

Cat. No.: GA23380

### Chemical Properties

Cas. No.

Formula  $C_{44}H_{69}N_{11}O_{13}$  M.Wt 960.09

Solubility Soluble in DMSO Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Riteau et al. carried out immunofluorescence analysis using a rabbit antiserum directed against the nonapeptide VGGQDKPFL corresponding to the N-terminal sequence of KIR2DL4. The peptide did not discriminate between p49 and KIR2DL4 receptor.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA