

---

**Product Data Sheet**

---

Product Name: Proadrenomedullin (12-20) (human)

Cat. No.: GA23385

**Chemical Properties**

Cas. No. 186027-43-0

Formula  $C_{56}H_{86}N_{18}O_{11}$  M.Wt 1187.41

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The proadrenomedullin fragment 12-20 has vasodepressor activity in the rat and cat. Proadrenomedullin (12-20) (PAMP 12-20) was approximately 3-fold less potent than the full-sequence peptide, but the duration of the vasodepressor response to this fragm

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA