

---

## Product Data Sheet

---

Product Name: Prokineticin 2 Isoform 2 (human)

Cat. No.: GA23388

### Chemical Properties

Cas. No. 423206-00-2

Formula  $C_{379}H_{606}N_{114}O_{101}S_{13}$

M.Wt 8792.55

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Prokineticin 2 (Bv8) is one of the first neuropeptides that appears to directly convey circadian rhythm signaling to other CNS structures and, therefore, it may be a critical output molecule regulating biological rhythms. The peptide has the capacity t

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA