

---

**Product Data Sheet**

---

Product Name: pTH-Related Protein Splice Isoform 3 (140-173) (human)

Cat. No.: GA23418

**Chemical Properties**

Cas. No. 139872-85-8

Formula  $C_{186}H_{313}N_{53}O_{44}S_2$ 

M.Wt 4059.99

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA