

---

**Product Data Sheet**

---

Product Name: rec IGF-I (human)

Cat. No.: GA23454

**Chemical Properties**

Cas. No. 67763-96-6

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

approx. ED = 0.8 ng/mL IGF-I regulates the proliferation and differentiation of a multitude of cell types, and also exhibits insulin-like effects, which are, in the case of glucose metabolism, opposite to those of the growth hormone (GH). IGF 'replacem

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA