

---

**Product Data Sheet**

---

Product Name: rec Leptin (human)

Cat. No.: GA23464

**Chemical Properties**

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Leptin, an adipose tissue-derived protein encoded by the obese gene, has been implicated in the control of food intake, energy expenditure and reproduction. Additional roles in the regulation of immunity, inflammation, and hematopoiesis have been proposed.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA