
Product Data Sheet

Product Name: rec TNF- α (human)

Cat. No.: GA23474

Chemical Properties

Cas. No. 94948-59-1

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This monocyte-macrophage-derived cytokine was first characterized by its in vivo antitumor properties and in vitro cytotoxicity to certain transformed cell lines. TNF- α also exerts a growth-enhancing effect on various cells and thus is a bifunctio

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA