

---

**Product Data Sheet**

---

Product Name: Relaxin H2 (human)

Cat. No.: GA23477

**Chemical Properties**

Cas. No. 99489-94-8

Formula C<sub>256</sub>H<sub>408</sub>N<sub>74</sub>O<sub>74</sub>S<sub>8</sub>

M.Wt 5963.04

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Relaxin H2 is the major stored and circulating isoform of the peptide hormone relaxin in humans. This pregnancy hormone has also shown cardioprotective activity. The observed long-term effects of relaxin H2 on connective tissues may be due to alteratio

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA